



Multi-Pan-Coater

GMPC



We set the standards

Glatt Multi-Pan-Coater

The new Glatt Multi-Pan-Coater (GMPC) is the newest development of the traditional Glatt Drum Coater. The GMPC is designed specifically for research and development for coating tablets with aqueous or organic solvent solutions.

As an option, it is also possible to coat cores with sugar syrup. The Multi-Pan has a compact construction and is suitable for GMP, features making it especially appropriate as a universal Coater.



All necessary components including the air intake and exhaust air system, are integrated into the machine housing. The Multi-Pan can be easily installed as a pilot machine.

The PLC control system records all necessary parameters, enabling the system to be scaled up to production plant size without any difficulty.




Due to the ease of the drum change, batch sizes from 2.0 to 56,0 litres can be processed.



Technical datas

MULTI-PAN type GMPC II

Technical data on request



Addresses

**The Glatt company
employs over
1500 people in
20 countries.**



Glatt® Maschinen- & Apparatebau AG

Kraftwerkstrasse 6
4133 Pratteln 1/Switzerland
Phone: +41 61 8 26 47 47
Fax: +41 61 8 26 48 48
E-mail: glatt-ag@glatt-ag.ch

Glatt® INOX AG

Wöschnauerstrasse 38
5012 Schönenwerd/Switzerland
Phone: +41 62 8 58 40 00
Fax: +41 62 8 58 40 10
E-mail: salesoffice@inox-glatt.ch



Glatt® GmbH Process Technology

Bühlmühle
79589 Binzen/Germany
Phone: +49 7621 6 64 0
Fax: +49 7621 6 47 23
E-mail: mail@glatt.de



Glatt® Ingenieurtechnik GmbH

Nordstraße 12
99427 Weimar/Germany
Phone: +49 3643 47 0
Fax: +49 3643 47 12 31
E-mail: git@glatt-weimar.de



Glatt® Systemtechnik GmbH

Grunaer Weg 26
01277 Dresden/Germany
Phone: +49 351 25 84 0
Fax: +49 351 25 84 328
E-mail: info@glatt-dresden.de



Glatt® Norden ApS.

Hejrevang 24
3450 Allerød/Denmark
Phone: +45 48 14 22 44
Fax: +45 48 14 22 55
E-mail: info@glattnorden.dk



Glatt® Pharmatech S.à.r.l.

Parc Technologique
6, rue Louis Néel
21000 Dijon/France
Phone: +33 3 80 74 32 64
Fax: +33 3 80 74 43 95
E-mail: glatt-pharmatech@glatt-pharmatech.fr



Glatt® Air Techniques Inc.

20 Spear Road
Ramsey, NJ 07446 USA
Phone: +1 201 8 25 87 00
Fax: +1 201 8 25 03 89
E-mail: info@glattair.com



Glatt® Protec Ltd.

Swannington Rd. Cottage Lane Ind. Est.
Broughton Astley
Leicester LE9 6TU/Great Britain
Phone: +44 1455 28 58 58
Fax: +44 1455 28 55 10
E-mail: info@glatt-protech.co.uk

Glatt® Labortecnic SA

C/ Cuzco, 55
08030 Barcelona/Spain
Phone: +34 93 3 45 67 00
Fax: +34 93 3 46 13 04
E-mail: glatt_labortec@cambrabcn.es

www.glatt.de



**We set the
standards**