

Fluid Bed Systems



MINI-GLATT



We set the standards

Mini-Glatt

Glatt, the global market leader in fluid bed technology, is pleased to present a very compact laboratory unit.

The mini fluid bed machine "Mini-Glatt" processes minimum quantities with maximum quality and is therefore ideal for drying or granulation trials on a very small scale.

This innovative small machine is easy to handle and extremely easy to clean. These properties make it perfect for use in pharmaceutical research and development laboratories.

The small and yet technically refined mini fluid bed machine offers convincing features: Small footprint, integrated and turnkey technology.

The special flow bed guarantees perfect fluidization of even extremely difficult products in extra-small batches.

Glatt developed a new spray nozzle specially for spraying minimum fluid volumes. Patented metal filter cartridges are used as outlet air filters.



The "Mini-Glatt": ideal for drying or granulation trials on a very small scale



The attractive price of the "Mini-Glatt" system won't strain budgets thanks to its low investment cost, while the small product application volumes help to save costs.



Technical data

Compliance
regulation
(TEX 100a)

Technical data on request

The Mini-Glatt is extremely easy to handle: without any difficulty, it can be stripped down and cleaned in seconds. All parts which come in contact with the product are made of high-grade steel (AISI-316L).

For proces-
ing solvents



**Trouble-free
handling**



Filter housing



Metal filter cartridges



Spray/Transition housing



Product bowl with lower plenum



Wurster insert

The Mini-Glatt can be stripped down in seconds with very little effort



Addresses

**The Glatt company
employs over
1,500 people in
20 countries.**



Glatt® GmbH Process Technology

Bühlmühle
79589 Binzen/Germany
Phone: +49 7621 6 64 0
Fax: +49 7621 6 47 23
E-mail: mail@glatt.de



Glatt® Ingenieurtechnik GmbH

Nordstraße 12
99427 Weimar/Germany
Phone: +49 3643 47 0
Fax: +49 3643 47 12 31
E-mail: git@glatt-weimar.de



Glatt® Systemtechnik GmbH

Grunaer Weg 26
01277 Dresden/Germany
Phone: +49 351 25 84 0
Fax: +49 351 25 84 328
E-mail: info@glatt-dresden.de



Glatt® Maschinen- & Apparatebau AG

Kraftwerkstraße 6
CH-4133 Pratteln/Switzerland
Phone: +41 61 8 26 47 47
Fax: +41 61 8 26 48 48
E-mail: glatt-ag@glatt-ag.ch



Glatt® Air Techniques Inc.

20 Spear Road
Ramsey, NJ 07446 USA
Phone: +1 201 8 25 87 00
Fax: +1 201 8 25 03 89
E-mail: info@glattair.com



Glatt® Pharmatech S.à.r.l.

Parc Technologique
6, rue Louis Néel
21000 Dijon/France
Phone: +33 3 80 74 32 64
Fax: +33 3 80 74 43 95
E-mail: glatt-pharmatech@glatt-pharmatech.fr

Glatt® Norden ApS.

Hejrevang 24
3450 Allerød/Denmark
Phone: +45 48 14 22 44
Fax: +45 48 14 22 55
E-mail: info@glattnorden.dk

Glatt® Protech Ltd.

Swannington Rd. Cottage Lane Ind. Est.
Broughton Astley
Leicester LE9 6TU/Great Britain
Phone: +44 1455 28 58 58
Fax: +44 1455 28 55 10
E-mail: info@glatt-protech.co.uk

www.glatt.de



**We set
the standards**